



RBP Chemical Technology, Inc.
150 S. 118th Street • P.O. Box 14069
Milwaukee, WI 53214-0069
www.rbpchemical.com
Corporate: 414-258-0911
Fax: 414-258-7908
RBP East: 508-865-0800

Product Data Sheet

Date: 08/01/11
Supersedes: 09/07/10
PRODUCT #: J5049

AQUA SERIES 180

Two Step Acidic Sheetfed Fountain Concentrate

Description:

A two step low VOC acidic fountain concentrate designed to eliminate calcium build-up on the rollers and keep calcium from depositing on blankets and plates. **AQUA SERIES 180** contains desensitizers and wetting agents that provide excellent ink and water balance on all types of sheetfed presses.

Directions for Use:

Mix 3 oz/gallon of **AQUA SERIES 180** to each gallon of water. Increase dosage if plate sensitivity is occurring and decrease dosage if ink emulsification is a problem. Each oz/gallon of **AQUA SERIES 180** will increase conductivity by 700 micromhos/cm. Use **AQUA SERIES 180** with an RBP alcohol replacer such as **AQUANOL® 600**, to eliminate the use of alcohol.

Specifications:

Appearance:	Hazy blue green liquid
Density:	9.0 lbs/gal
Conductivity:	700 micromhos/cm per ounce
pH, diluted:	3.5-4.5
VOC Content:	0.7 lbs/gal, 84 g/L, 7.8 %

This product should be used only for its intended purpose. The information stated above is based on our laboratory tests and experience, and is accurate to the best of our knowledge. Since actual use is beyond our control, the recommendations or suggestions are made without warranty, expressed or implied.

